

Helicobacter pylori Stool Antigen Notice 幽門螺旋菌大便測試須知

Preparation and Instruction for Investigation 測試準備及須知:

- 1. Patient should have bismuth salts or proton pump inhibitors stopped 14 days prior to the testing. 病人需於 14 日前停用鉍製劑、氫離子幫浦阻斷劑
 - a. Bismuth preparation (Pepto Bismol) 鉍製劑 (鹼式水楊酸铋)
 - b. Proton pump inhibitors (e.g.: Losec, Prevacid, Pantoloc, Nexium, Pariet). 氫離子幫浦阻斷劑 (例如: 樂酸克、蘭索拉唑、保衛康治潰樂腸膜溶衣錠、耐適恩錠、百抑潰) 。

For post eradication evaluation, the test should not be carried out earlier than four weeks after active treatment has been stopped.

對於治療後的病人則需在療程結束 4 周後才進行測試。

- 2. Do not use any medication e.g. laxatives, to help pass the stool. 不要使用任何藥物,如瀉藥來幫助排便。
- 3. Urinate into the toilet if needed. 如有需要,可先小便至馬桶內。
- 4. Place sheets of plastic wrap over the toilet bowl or use a clean bowl or bedpan. Do not let urine or water touch the stool specimen. 將塑膠膜撲在廁缸上,中心處留出低窪處或使用乾淨的廁缸或便盆,請勿讓尿液或水接觸糞便樣本。
- 5. Use the spoon attached to the cap to take a portion of stool into a clean container without preservative. 使用連在蓋子上的小勺,取適量糞便放入塑膠容器中。不可使用盛放清潔劑或防腐劑的容器。
- 6. Wash hands with soap and water. Write your Full Name and Time of collection on container, place the container into plastic bag and seal. 用肥皂和清水洗手。在容器上註明姓名及採集日期和時間,將容器放入塑膠袋中並密封。
- 7. Bring the container back to the Clinical Laboratory within 2 hours of collection. If a delay of more than 2 hours is expected, please refrigerate the sample. 在採集2小時內將容器帶回化驗室,延遲可能會影響檢測結果。如樣本未能在採集後2小時內送回 化驗室,請把標本冷藏。

Pathology Department Version 4.0

Effective Date: 13 May 2024